

Product Name :
Educational Lab Glass Rod

Product Code :
NLAB-EDUSCHOOLEQ340107



Description :

Educational Lab Glass Rod

Technical Specification :

Solid glass rod dia 5mm x 170mm.

Laboratory glass rods are used in various applications like mixing of chemicals, heating experiments, magnetic charge conductions etc.

Sturdy glass material is used in the manufacturing of rods.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003