

Product Name :
Educational Lab Haemoglobin Meter, Microprocessor Based

Product Code :
NLAB-EDUSCHOOLEQ340137



Description :

Educational Lab Haemoglobin Meter, Microprocessor Based

Technical Specification :

Microprocessor Haemoglobin Meter is an indispensable need of medical laboratories.

This meter is used to measure the haemoglobin concentration in a blood sample.

The measurement haemoglobin concentration is carried out at wavelength of 546 nm.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003