

Product Name :
Educational Lab Insulated Copper Wire

Product Code :
NLAB-EDUSCHOOLEQ34052



Description :

Educational Lab Insulated Copper Wire

Technical Specification :

Copper strips find many applications as electrical conductors, switch terminals, hardware used in electrical circuits, electrical conductors.

These Insulated Copper Wire are available in hard, half-hard and soft conditions.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003