

**Product Name :**  
Educational Lab Insulated Copper Wire

**Product Code :**  
NLAB-EDUSCHOOLEQ34052



**Description :**

Educational Lab Insulated Copper Wire

**Technical Specification :**

Copper strips find many applications as electrical conductors, switch terminals, hardware used in electrical circuits, electrical conductors.

These Insulated Copper Wire are available in hard, half-hard and soft conditions.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003