

Product Name :
Educational Lab Kipps Apparatus

Product Code :
NLAB-EDUSCHOOLEQ340128



Description :

Educational Lab Kipps Apparatus

Technical Specification :

This apparatus is made of polyethylene has three main portions.

The top portions is fitted with a thistle funnel and a bulb structure through which dilute acid can be poured will reacts with a metal sulphide in the base of the Apparatus.

HS gas then evolved collects in the middle dome from where it can be released through a leakproof stopcock fitted with a self lubricating PTFE plug.

The base is provided with a residual outlet fitted with a cap for easy cleaning.

Particulars:

250mm

500mm

1000mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003