


Product Name :
Educational Lab Lid for Copper Calorimeter

Product Code :
NLAB-EDUSCHOOLEQ34056



Description :

Educational Lab Lid for Copper Calorimeter

Technical Specification :

Copper calorimeter with bakelite lid at the top fitted with rubber stopper for thermometer and metal stirrer.

Size: 100 x 75mm (Height x Diameter).

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003