

Product Name :
Educational Lab Liebig Condenser

Product Code :
NLAB-EDUSCHOOLEQ340130



Description :

Educational Lab Liebig Condenser

Technical Specification :

Liebig condenser is manufactured to condense vapours through a centre tube.

It has 24/40 joints, and 300 mm jacket length.

The material in which it is made is borosilicate glass and is shaped through hand blowing technique.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003