

Product Name :
Educational Lab Magnascope

Product Code :
NLAB-EDUSCHOOLEQ340132



Description :

Educational Lab Magnascope

Technical Specification :

Lens : Dia 150 mm made from fine optical glass

Effective Lens Dia. : 135mm

Magnification : 3.5x

Body Material : Aluminum & C.I.

Body Color : Siemens Grey.

Top Illumination : 22W Circular Fluorescent Tube Dia.216mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

