

Product Name :
Educational Lab Micropipettes Multi Channel

Product Code :
NLAB-EDUSCHOOLEQ340144



Description :

Educational Lab Micropipettes Multi Channel

Technical Specification :

Multi-channel pipettes with smart soft finger-grip of Thermoplastic Elastomer ensure minimum fatigue of hand and fingers.

It has a soft force tip loading and ejection.

Manifold can rotate 360° for easy right or left hand operation.

It also allows consistent sample loading.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003