


**Product Name :**  
Educational Lab Vernier Caliper

**Product Code :**  
NLAB-EDUSCHOOLEQ34029



**Description :**

Educational Lab Vernier Caliper

**Technical Specification :**

Optional graduations: 0.02mm,0.02mm/0.001", 0.05mm.

With thumb knob.

monoblock, made of hardened stainless steel, satin chrome plated reading surface.

Round depth rod is optional.

Titanium plated beam is optional.

Packed in leathery pouch or plastic box.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003

