

Product Name :
Elastic Recovery Mould

Product Code :
NLAB-BTMENLAB150015



Description :

Elastic Recovery Mould

Technical Specification :

A set of three moulds as per IRC-SP-53-2002 made of brass /gunmetal with the shape, dimensions and tolerances. The ends B & B` are known as clips, and the parts A & A' are known as sides of the mould. The dimensions of the mould shall be such that when properly assembled it will form a briquette specimen having the following dimensions. Total length : 75 ± 0.5 mm

A = 36.5 ± 0.1 MM

B = 30.0 ± 0.1 MM

C = 17.0 ± 0.1 MM

D = 10.0 ± 0.1 MM

E = 10.0 ± 0.1 MM

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments