

Product Name :
Electric Backrest

Product Code :
NLABMEDICAL-M2560002



Description :

Electric Backrest

Technical Specification :

Specification :-

- * Frame is made of Mild Steel Tube.
- * Top uniformly perforated CRCA sheet.
- * Back rest adjustable on Electrically Operated with Remote.
- * Supply with foam mattress 15 cm thick.
- * Size: 89 L x 74 W cm.
- * Finish in epoxy powder coated.

We are Manufacturer, Exporters and Suppliers of Electric Backrest at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003