

Product Name :
Electric C-Arm Table with Top Slide

Product Code :
NLABMEDICAL-M3080001



Description :

Electric C-Arm Table with Top Slide

Technical Specification :

* This Electric sliding top O.T Table is suitable for all kind of surgeries, i.e. general surgery, vascular surgery, cardio, neurology, urology, gynecology, proctology, laparoscopy, traumatic surgery, plastic surgery and others.

* The base & column of the table is made of high quality medical stainless steel (Grade – 304).

* Non- reflecting surface, antibacterial & easy to clean.

* We provide with a top sliding of 5inch back/forward.

* Height adjustment, Lateral Tilt, Trendelenburg / Reverse Trendelenburg and Flex – Reflex, Chair Position and Longitudinal Sliding top are realized by Remote control.

Technical Specs:-

Length of Table 1900 mm

Width of Table Top 520 mm

Minimum Height (without Mattress) Min. 750 mm

Maximum Height (without Mattress) Max. 1000 mm

Trendelenburg 25°

Reverse Trendelenburg 25°

Lateral Tilt $\pm 20^\circ$

Back Section $\pm 70^\circ$

Head Section $\pm 60^\circ$

Leg Section 0°-90°

Top Sliding 5° Reverse & 5° Forward

Floor Locking System Center Floor Locking

We are Manufacturer, Exporters and Suppliers of Electric C-Arm Table with Top Slide at tender and export competitive prices from India

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003