

Product Name :
Electric Examination Couch

Product Code :
NLABMEDICAL-M2540003



Description :

Electric Examination Couch

Technical Specification :

Specification :-

Frame Work of the Couch is made of Rectangular & square tube.

It has 2 Section Cushioned top.

The body section is mounted on 12.5 cm castors.

The Back Section is adjustable by electric actuator control by remote.

The height is adjustable by electric actuator control by remote upto 80cm.

Product Size : 183 L x 60 W x 60-80 H cms.

Knock Down Construction to save freight.

We are Manufacturer, Exporters and Suppliers of Electric Examination Couch at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003