

Product Name :
Electric Obstetric Table

Product Code :
NLABMEDICAL-M2540001



Description :

Electric Obstetric Table

Technical Specification :

Obstetrics operating room table perfect for all kinds of gynecological operations.

The cover for base and column is made of 304 stainless steel, good performance on anti-rust and easy to clean.

Technical Data

Length of Table : 1600 mm.

Width of Table Top : 550 mm.

Minimum Height (without Mattress) : Min. 750 mm.

Maximum Height (without Mattress) : Max. 1000 mm.

Lateral Tilt : $\pm 20^\circ$.

Back Section : $\pm 70^\circ$.

Floor Locking System : Manual.

Features

Functional and versatile.

Flexibility and precision.

Hygienic design.

Easy cleaning.

User- friendly.

STEEL QUALITY 304 MEDICAL GRADE STAINLESS STEEL.

Standard Accessories

- * 'L' Shaped Anaesthetic Screen : 1 Pc.
- * Padded Shoulder Support : 1 Pair.
- * Radio-translucent Armboard with pad : 2 Pcs.
- * Goepel Knee Crutch with straps : 1 Pair.
- * PU Foam Mattress : 1 Set.
- * Water proof table top cover : 1 Pc.
- * S.S. Douche Tray with drain pipe : 1 Pc.

We are Manufacturer, Exporters and Suppliers of Electric Obstetric Table at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003