

**Product Name :**  
Electrical Resistance Tester Machines

**Product Code :**  
NLAB-TESTINGLABM37001



**Description :**

Electrical Resistance Tester Machines

**Technical Specification :**

- Reliable functioning
- Smooth operation
- Precisely engineered

**Specifications:**

- Test voltage :  $100 \pm 2$  volts DC
- Internal electrodes : 5 mm diameter steel balls (total mass 4 kg)
- External electrode : 100 mm x 230 mm x 6 mm thick copper plate
- Voltmeter : 199.9 volts DC (digital)

Ammeter : 199.9 mA, 1.999 mA, and 19.99 ?A (digital)

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003