

Product Name :
Electrical test screwdriver

Product Code :
NLAB-ELECTRICALAB210005



Description :

Electrical test screwdriver

Technical Specification :

Specifications:

- Blade made of Cr-V. 100-500V AC.
- With clip.
- Super-long life Neon bulb can be used over 20,000 times.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003