

**Product Name :**  
Electro Hydraulic Trainer

**Product Code :**  
NLAB-ENGINEERINGLB101002



**Description :**

Electro Hydraulic Trainer

**Technical Specification :**

Electro Hydraulic Trainer

It is precisely engineered using premium quality raw material and modern machinery at our unit in adherence with the prevailing industrial norms and standards.

Our experts conduct various quality inspections to ensure that Electro Hydraulic Trainer stands high on strength, quality, effectiveness and has a long life.

Technical Specification:-

Operating pressure 40 bar approx., Flow rate 4-8 lpm approx.

Hydraulic Tank capacity 20 liters.

Geared pump coupled with induction motor 2 HP 220 V AC,50 Hz, and mounted externally over the tank with anti vibration pads.

4/3 Double solenoid valve AC operated.

4/2 Single solenoid valve AC operated.

4/3 Way hand lever valve (Manual).

Double acting cylinder (Stroke length 100 mm & diameter 32 mm).

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003