

Product Name :
Electrolysis Apparatus

Product Code :
NLAB-EDUSCHOOLEQ143002



Description :

Electrolysis Apparatus

Technical Specification :

A simple type of Electrolysis Apparatus, but demonstrates the principle as efficiently as other expensive apparatus.

Consists of a plastic container 3 inch dia x 2.5 inch ht, total height 3.75 inch, with mounted carbon electrodes.

The bottom of the container has cut outs so that leads can be connected to the electrodes with ease.

Supplied complete with two alligator clip leads (1 red, 1 Black), two graduated test tubes and instruction manual.

Works on 12v dc.

Measurements are approximate.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003