

Product Name :
Electromagnetic Digital Sieve Shaker

Product Code :
NLAB-SOILELAB140040



Description :

Electromagnetic Digital Sieve Shaker

Technical Specification :

Electromagnetic Digital Sieve Shaker performs a vertical sieving motion and by removing the upper plate it can also be used as a vibrating table for one concrete specimen. The vertical sieving motion is provided by a very effective electromagnetic unit which is designed to obtain the best results with sand and aggregates. Supplied complete with timer.

Sieve Capacity

12 pieces of 200 mm (8") sieves + pan and cover

8 pieces of 300 mm (12") sieves + pan and cover

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003