

Product Name :
Electromatic O.T. Table

Product Code :
NLABMEDICAL-M2760001

Description :

Electromatic O.T. Table

Technical Specification :

It is equipped with imported electric actuators for smooth & effortless working of the table.

It ensures the actuators working in a constant speed without any noise & is designed to run even at low voltage to carry on the surgery.

Key Structure:

The table has been made from acid proof, stainless steel (the face of the table is easy to clean and immune to disinfection agents)

Strong & rust free base of stainless steel is provided to make the working conditions pleasant for the doctors & nurses.

Orthopedic attachment (hips and legs support) including the patient's pelvis support and legs support enables the limbs strain performing.

Design of the Orthopedic attachment & eccentric column provides monitoring of the limbs with the C-ARM during the operation.

The various positions such as Up-Down movements, Trendelenburg Reverse Trendelenburg.

Lateral Tilts, Flex- Reflex, Chair Position as realized by remote control.

Leg section & head section are interchangeable.

Technical Constraints:

Top Dimension: L 1930 XW 553 mm

Height Adjustment: 762 mm -1012 mm

Trandelenburg / Reverse : 30° / 25°

Lateral Tilt: 20° / 20°

Kidney Elevator: 150 mm

Back Rest (up / Down): 80° / 25°

Leg Rest (up / Down): 15° / 90°

Head Rest (up / Down): 20° / 60°

Power Supply: 24v Dc

Battery Backup (Optional): 2-3 Hours

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003