

**Product Name :**  
Electronic Circuit Kit

**Product Code :**  
NLAB-SCHOOLABS16001



**Description :**

Electronic Circuit Kit

**Technical Specification :**

Speed LED

flashing light, EK133 6-LED disco

light,3-step flashing light, 19 LED starburst,

nurse siren (European sound) ready. Loudspeaker (Hospital siren)

Musical door bell with speaker

Wake-up with speaker

Bell sounds with buzzer (Bell sounds)

Digital laughter with speaker

Christmas music with speaker

siren 6 Sound with loudspeaker

Dimmable 500 w

Fan dimmed 1000 w

Water pump alarm and cut-off circuit (full, empty)

intruder detection with loudspeaker.

electronic guard with speaker

AM radio IC system with speaker

Simple AM â€œradio with earpiece

3-level water

indicator Rain alarm with speaker (9V)

overflow alarm with speaker

refrigerator door closed warning EK912

pressure meter Light intensity 5 levels

mosquito repellent

cycle repellent mice – cockroaches

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003