

Product Name :
Electronic Leak Detector

Product Code :
NLAB-VOCATIONALTVET13016



Description :

Electronic Leak Detector

Technical Specification :

Operating Temp: 32° to 125°F

Power Supply: 2 C-Size 1.5V Ni-MH Battery

Battery Life: 8 hrs per cycle (1000 cycles)

Warm-Up Time: 25 Seconds

Ultimate Sensitivity: 1/10 oz per year

Probe/Core Length: 15"

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

