

Product Name :
Electronic Stimulator

Product Code :
NLAB-SCIENTIFICEQ28009



Description :

Electronic Stimulator

Technical Specification :

Available in four models:

Student Stimulators Transistorised : Out put 0-25 DC in two steps, duration .5 to 1.5ms frequency 5 to 100 c/s in 11 steps, modes, single, repetitive or an external trigger, for 220V 50Hz.

Student Electronic Stimulator Latest Design :

Stimulator General Purpose : Provides rectangular pulse of controlled amplitude 0.1 – 100V. Frequency 0.1 – 1000 PPS & width 0.01 to 100mm.

Research Stimulator With R.F. Isolation unit : Provides rectangular pulse of controlled polarity, Amplitude 1.5mV to 150V, frequency 0.1 to 10 KCS delay 10 micro sec to 1 sec & width 10 micro sec to 1 sec.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

