

Product Name :
Electronic Thermometer

Product Code :
NLAB-LBEMALAB230027



Description :

Electronic Thermometer

Technical Specification :

This series of thermocouple thermometers are micro controller based testers. With very intuitive interface, the user can fully operate the tester without second thought. They are digitally calibrated for long time stability. The whole series are equipped with PC interface for advance applications.

Features:

Timer function

PC interface

Auto power off

Resolution 0.1°C; 0.1°F

Backlit LCD

4 input channels

4 readings display

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003