

**Product Name :**  
Electronic Universal Testing Machines

**Product Code :**  
UMT-0003



**Description :**

Electronic Universal Testing Machines

**Technical Specification :**

**Features :**

- Motorised Control Valve.
- Load accuracy as high + 1%.
- Straining at variable speeds to suit a wide range of materials.
- Printer & PC graphs enable study the behavior of the material.
- Motor driven threaded columns for quick effortless adjustment of middle crosshead to facilitate rapid fixing of test specimen.
- Simplicity in reading because of digital readouts.
- Wide range of standard and special accessories.

- Easy change from plain to threaded and screwed specimens.
- Large effective clearance between columns enable testing of standard specimens as well structures.
- Simple controls for ease of operation.
- Robust straining frame of an extremely rigid construction.
- Fully enclosed and protected pressure transducer.
- RS 232 serial port for PC control of m/c.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-0171-2643080, 2601773