

**Product Name :**  
Electrophoresis Apparatus

**Product Code :**  
NLAB-EDUSCHOOLEQ45033



**Description :**

Electrophoresis Apparatus

**Technical Specification :**

The horizontal electrophoresis apparatus was designed with safety and convenience in mind.

Each casting tray with embossed gel ruler has rubber end caps and multiple slots for easy and versatile comb placement.

The injection-molded UV transparent chamber is made of thick acrylic, making it durable and leak-proof while the lid fits securely in place and can easily be removed.

It allows two 7x7cm gel trays to run simultaneously on one apparatus.

This unit requires 500mL of buffer and a power supply.

Electrophoresis units feature adjustable leveling feet, recessed electrical wires, and a safety stop that prevents the gel from running when the cover is not securely fitted.

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com  
**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

## Naugralabequipments