


Product Name :
Electrophoresis Classroom Kit

Product Code :
NLAB-EDUSCHOOLEQ45034



Description :

Electrophoresis Classroom Kit

Technical Specification :

The injection-molded UV transparent chamber is made of thick acrylic, making it durable and leak-proof while the lid fits securely in place and can easily be removed.

Each casting tray with embossed gel ruler has rubber end caps and multiple slots for easy and versatile comb placement.

Electrophoresis units feature adjustable leveling feet, recessed electrical wires, and a safety stop that prevents the gel from running when the cover is not securely fitted.

The horizontal electrophoresis apparatus was designed with safety and convenience in mind.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003