

Product Name :
Electrophoresis Dual Station

Product Code :
NLAB-EDUSCHOOLEQ45032



Description :

Electrophoresis Dual Station

Technical Specification :

Electrophoresis units feature adjustable leveling feet, recessed electrical wires, and a safety stop that prevents the gel from running when the cover is not securely fitted.

The horizontal electrophoresis apparatus was designed with safety and convenience in mind.

The injection-molded UV transparent chamber is made of thick acrylic, making it durable and leak-proof while the lid fits securely in place and can easily be removed.

Supports up to 16 students.

Each casting tray with embossed gel ruler has rubber end caps and multiple slots for easy and versatile comb placement.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

