

Product Name :
Electrostatic Kit

Product Code :
NLAB-VOCATIONALTJET30016



Description :

Electrostatic Kit

Technical Specification :

1 each reel nylon.

Cellulose acetate strip.

Polythene strip.

Wire stirrup.

Rubber cloth.

Wire hook.

Electrophorus on handle and proof plane.

2 each of polythene tiles and aluminum cans.

We are well-known manufacturers, OEM suppliers of Electrostatic Kit for Physics Workshop Labs. Contact us for high quality Electrostatic Kit for Physics Workshop Labs for schools lab, college lab, universities, research labs,

various teaching and workshop training laboratories and industries in India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003