

<b>Product Name :</b> Emergency and Recovery Trolley (Hydraulically Operated)	<b>Product Code :</b> NLABMEDICAL-M2770001
<b>Description :</b>	
Emergency and Recovery Trolley (Hydraulically Operated)	
<b>Technical Specification :</b>	
Framework made of round & rectangular tubes mounted on heavy-duty castors of 12.5 cm dia with two brakes.	
Offers trend enberg & reverse position and height adjustable by hydraulic pump fitted on one side of the trolley.	
Removable type Radio-lucent top in two sections for obtaining backrest and fitted with a cassette holding tray underneath.	
Supplied with swing away type safety railing telescopic saline rod adjustable, steel storage tray and oxygen cylinder cage.	
Finish in epoxy powder coated.	
Size: 210 L x 68 W x 65-85 H	

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003