

Product Name :
Encoder And Decoder Trainer

Product Code :
NLAB-ELECTRONICSAB170003



Description :

Encoder And Decoder Trainer

Technical Specification :

Encoder And Decoder Trainer For Electrical Lab Training.

- Study of Encoder & Decoder
- input switches with LED display
- Logic status indicator
- Set of Patch chords
- User's Manual
- House in an attractive ABS Plastic cabinet with cover

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003