

Product Name :
Engine Performance Trainer

Product Code :
NLAB-ENGINEERINGLB39001



Description :

Engine Performance Trainer

Technical Specification :

Engine Performance Trainer Lab Equipments Manufacturer, Suppliers & Exporters

Engine Performance Trainer

Features of Engine Performance Trainer

1. The Engine Performance Trainer is constructed from a late model “recycled” vehicle (FrontWheel-Drive or Rear-Wheel-Drive)

2. Includes all of the Following Automobile Systems for Study and Diagnosis:

Engine (4 to 8 Cylinder, Gasoline or Diesel)

Transmission (Automatic or Manual)

Heat/Ventilation/Air Conditioning

Electrical

Exhaust

Engine Cooling

Emissions

3. Aids in Developing Real-World System Skills and Knowledge

4. Performs Actual Service Manual Diagnostic Procedures

5. Includes an ATech Fault/Simulator Panel with Multiple Faulting Options and Features Including:

Real-Time voltage readouts for Various Selectable Sensors

Individual or Multiple Sensors are selectable for Voltage Substitution

Voltage Adjustment Potentiometers for each Sensor

Interactive Demonstration/Study of Fuel Control by Oxygen Sensor Feed back

Demonstrates Catalytic Convertor "Punch-Through"

Investigates Downstream Oxygen Sensor's Control of Fuel Trim

6. Breakout Box with Test Points for every PCM Pin

7. Data Link Connector (DLC) for Scan Tool Diagnostics

8. Trainer Mounted on Casters for Ease of Mobility

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003