

Product Name :
ENT Patient Chair

Product Code :
NLABMEDICAL-M3220019

Description :

ENT Patient Chair India, ENT Patient Chair manufacturer India, suppliers, Medical ENT Patient Chair India, Hospital ENT Patient Chair India

Technical Specification :

Fully Flexible, Rotating Recliner Chair

Noiseless, DC-Motor driven, Up & Down movements and inclination

Board-based design for high stability

Comfortable polyurethane-filled seat

Arm rest with swing-up option

Adjustable footrest, which can be fully raised to Sleeping position

Under-calf cushioning in sleeping position

Ideal movement control from chair side top control panel

Foot Switch for fluid up / down movement control

Optimal overall ergonomics

Minimum measurement needed ~8/10 sq. ft.

Power Supply-220V/50HZ/ 5 Amp

Seat Minimum Height- 550mm

Seat Maximum Height- 710mm

Up & Down Stroke- 160mm

Seat Rotation- 0 to 180 deg.

Back & Forward Movement- 95° to 750°

Power Consumption- 600mA

Net Weight- 200 Kgs.

Load Lifting- 200 Kgs.

Motor Available- 230V AC-24V Dc

Minimum are required- 8 feet to 10 feet

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003