

Product Name :
Equipment For Studying Rotational Motion India

Product Code :
NLAB-VOCATIONALTVET30005



Description :

Equipment For Studying Rotational Motion

Technical Specification :

Adapters

Pulley and rod

Chronometer

Springs

Turning platform, base support 22 x 22 cm and track 51 x 4 cm

Disk \varnothing 23 cm

Ring, outer \varnothing 12,5 cm, inner \varnothing 10,5 cm

2 trucks (1 mobile, 1 steady)

Masses from binary

DC motor, 24 V, 1500 rpm

Accessories for studying centrifugal force

Pulley from binary

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003