

Product Name :
Erichen Cupping Machine

Product Code :
NLAB-MECHANICALAB250009



Description :

Erichen Cupping Machine

Technical Specification :

We have gained expertise in manufacturing, supplying and trading high quality Erichsen Cupping Testing Machine. Offered range is designed to reveal the cupping qualities of metal sheets and strips. These technologically advanced machines are used for testing adhesion, elasticity and priority of coat of paints or varnish by way of comparison.

Features:

Low maintenance

Longer functional life

Excellent performance

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003