


Product Name :
Exhaust Gas Calorimeter

Product Code :
NLAB-MECHANICALAB210005



Description :

Exhaust Gas Calorimeter

Technical Specification :

Exhaust Gas Calorimeter: The experimental exhaust gas calorimeter water heat exchange to be used with multi-cylinder engine. Small vertical finned tubes enclosed in insulated stainless steel chamber. Exhaust gas from engine can be connected to the calorimeter and flow around the finned tubes while cooling water flows inside the tubes. Chamber can be disassembled for cleaning

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003