

Product Name :
Extension And Compression Of Springs Apparatus

Product Code :
NLAB-ENGINEERINGLB26003



Description :

Extension And Compression Of Springs Apparatus

Technical Specification :

Extension And Compression Of Springs Apparatus Lab Equipments Manufacturer, Suppliers & Exporters

Description of Extension And Compression Of Springs Apparatus:

Wall mounted apparatus (for loading spiral springs in either tension or compression with direct gravitational loads) with accurately divided scale and vernier. Supplied complete with four compression and four extension springs each approximately 18 cm long (7") long. Mass Hanger and Masses for (Hanger 1, Masses 13 x 1 kg)

Experimental Capability of Extension And Compression Of Springs Apparatus:

Determines load and deflection relationships for close coiled springs.

Modulus of rigidity for spring material.

Introduction to Spring Rate

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003