

**Product Name :**  
Extension And Compression Of Springs Apparatus

**Product Code :**  
NLAB-ENGINEERINGLB26003



**Description :**

Extension And Compression Of Springs Apparatus

**Technical Specification :**

Extension And Compression Of Springs Apparatus Lab Equipments Manufacturer, Suppliers & Exporters

Description of Extension And Compression Of Springs Apparatus:

Wall mounted apparatus (for loading spiral springs in either tension or compression with direct gravitational loads) with accurately divided scale and vernier. Supplied complete with four compression and four extension springs each approximately 18 cm long (7") long. Mass Hanger and Masses for (Hanger 1, Masses 13 x 1 kg)

Experimental Capability of Extension And Compression Of Springs Apparatus:

Determines load and deflection relationships for close coiled springs.

Modulus of rigidity for spring material.

Introduction to Spring Rate

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003