

**Product Name :**  
Fatigue Testing Machine

**Product Code :**  
TN366



**Description :**

Fatigue Testing Machine

**Technical Specification :**

Fatigue Testing Machine Light weight, compact size, simple design Table model Simple lever system for changing load 10 Digit electronic counter instead of mechanical counter can be supplied at additional Machine with bending movement up to 400kg/cm can be offered on request. Max Bending Moment: 400Kg. B.M adjustable:0 - 400 Kg. Specimen Dia.: up to 100mm Motor speed:2800 RPM

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003