

Product Name :
Faucet

Product Code :
NLAB-LABORATORYEM81004



Description :

Faucet

Technical Specification :

Application: Office Building, Hospital, School, Other

Material: Brass+PP

Function: 360 degree rotate

Usage: Control the effluent

360 degree rotation, convenient water discharge

Miniature spout,Dense water, not easy to splash, save water,

Coating: High-brightness and ultra-thick electroplating layer, corrosion-resistant, heat-resistant, and UV-resistant.

Valve core: adopts copper ceramic valve core, with a service life of 500,000 switching times, and a static pressure resistance of 35 bar

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003