

Product Name :
FEED Biasing of Transistor Trainer Digital Meters

Product Code :
NLAB-ELECTRICALAB210065



Description :

FEED Biasing of Transistor Trainer Digital Meters

Technical Specification :

SCOPE OF LEARNING:

Feed /Fixed Biasing of BJT

TECHNICAL SPECIFICATIONS:

Digital Meters:

Voltmeter 20VDC

Ammeter 200uA/2mA DC.

Power Supplies:

DC Supply IC Regulated +12V DC, 250mA.

Operated on Mains power 230V, 50Hz +10%

Components are mounted on the panels are:

Transistor

Required Resistors

SALIENT FEATURES:

Front panel built with high class insulated Printed Circuit Board sheet with well printed circuits and symbols.

Fuse for Short Circuit protection

Instruction manual.

Connections are brought out through 4mm Colored Sockets.

Patch Cords 4mm.

The trainer is housed in ABS Plastic cabinet.

Size of the trainer set 12"x8"

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003