

Product Name : Fetal Monitors	Product Code : NLABMEDICAL-M2740008
Description :	
Fetal Monitors	
Technical Specification :	
Highly sensitive.	
Wide 5.7 blue screen LCD.	
Full functional remote control.	
Self-define multi-event marks.	
Fetal movement auto detection.	
7-chip impulse ultrasound transducer.	
Built-in high-resolution thermal array printer.	
Auto power failure protection and reliable rechargeable.	
Preset usual functions shortcut keys for convenient operation.	
Signal intensity indication ensures transducer detecting position.	
High performance and self-adapting wide scope supply voltage design.	
Continuous monitoring of Fetal Heart Rate and Uterine Contraction (OP).	
User friendly interface with menu function settings, direct and ease-of-use.	

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com