


Product Name : Fetal Monitors	Product Code : NLABMEDICAL-M2740008
---	---




Description :

Fetal Monitors

Technical Specification :

Highly sensitive.
Wide 5.7 blue screen LCD.
Full functional remote control.
Self-define multi-event marks.
Fetal movement auto detection.
7-chip impulse ultrasound transducer.
Built-in high-resolution thermal array printer.
Auto power failure protection and reliable rechargeable.
Preset usual functions shortcut keys for convenient operation.
Signal intensity indication ensures transducer detecting position.
High performance and self-adapting wide scope supply voltage design.
Continuous monitoring of Fetal Heart Rate and Uterine Contraction (OP).
User friendly interface with menu function settings, direct and ease-of-use.



Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com