

Product Name : Fiber Optic Operating Light	Product Code : NLABMEDICAL-M3230009
Description : Fiber Optic Operating Light India, Fiber Optic Operating Light manufacturer India, suppliers, Medical Fiber Optic Operating Light India, Hospital Fiber Optic Operating Light India	
Technical Specification : Non-Drift light pipe with 360 degree rotation & lock on any position Variable Halogen light of High Intense 150w Bulb Idle for various Micro Surgery Operations Compact light head with Focussing Sleeve With one spare Bulb	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003