


Product Name : Fibre Optics Kit	Product Code : NLAB-VOCATIONALTVET30013
	
Description : Fibre Optics Kit	
Technical Specification : Data reception card Analog transmitter/receiver Data transmission card Receiver with loudspeaker Frequency generator Digital transmitter/receiver Transmitter with microphone	
Training Program:- Detecting signals, an analog receiver	

Transmitting a signal, a function generator

Transmission cation of sound, transformation of sound into electric signal, optical pulse, and their back conversion

Tyndall's experiment

Optical transmission between PCs (2 optional Personal Computers are required)

Examples of optical sensors (level sensor, force sensor, photogate sensor)

We are well-known manufacturers, OEM suppliers of Fibre Optics Kit for Physics Workshop Labs. Contact us for high quality Fibre Optics Kit for Physics Workshop Labs for schools lab, college lab, universities, research labs, various teaching and workshop training laboratories and industries in India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003