

**Product Name :**  
Field Monitoring Technology system

**Product Code :**  
NLAB-TESTINGLABM33001



**Description :**

Field Monitoring Technology system

**Technical Specification :**

Provides average wind speed, highest 3-second wind gust, and average wind direction for the measurement interval with a range of 0 to 170 mph (0 to 76 m/sec). A complete system requires an RX3000 Remote Monitoring Station, Wireless Manager, and Wireless Sensor.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003