

Product Name : Field Monitoring Technology system	Product Code : NLAB-TESTINGLABM33001
---	--



Description :
Field Monitoring Technology system
Technical Specification :
Provides average wind speed, highest 3-second wind gust, and average wind direction for the measurement interval with a range of 0 to 170 mph (0 to 76 m/sec). A complete system requires an RX3000 Remote Monitoring Station, Wireless Manager, and Wireless Sensor.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003