

Product Name :
Filing Cabinets 3 Drawer

Product Code :
NLABMEDICAL-M2570001



Description :

Filing Cabinets 3 Drawer

Technical Specification :

Filing cabinet is made of MS Sheet.

It has 3 Drawers.

It is fitted with Lock Key.

This is made of Mild Steel & Pre Treated Epoxy Powder Coated.

We are Manufacturer, Exporters and Suppliers of Filing Cabinets 3 Drawer at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003