

Product Name :
Film Strip-Slide Projector

Product Code :
NLAB-SCIENTIFICEQ15003



Description :

Film Strip-Slide Projector

Technical Specification :

Film Strip-Slide Projector is the latest, improved & versatile projector model. It is provided with efficient cooling system and heat resistant glass protection.

Key Points:-

- Moulded body
- Suitable to project 35mm film strips.
- 2? x 2? projection slides.
- 24V–150W projection illumination.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003